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
A photograph of a doctor in a light blue lab coat performing a laser procedure on a patient's face. The patient is lying down with their eyes closed. A CBS4 logo is visible in the bottom right corner of the image.

eCO₂ Laser

Used On:
Face, Neck, Chest

Can Improve:
Fine lines
Sun-damaged skin
Tone, texture

Print | Close this window

Lutronic(R) Becomes Official Laser for Lifestyle Lift(R)

Tue Jul 21, 2009 8:45am EDT

Lifestyle Lift(C) Selects eCO2(TM) as Its Laser of Choice

PRINCETON JUNCTION, N.J., July 21 /PRNewswire/ -- Lutronic, Inc., a leading developer of high quality, affordable medical and aesthetic laser systems, has signed a contract with Lifestyle Lift(C) to become their primary provider of laser therapy with its fractional CO2 laser system, eCO2(TM). Lifestyle Lift(C) is one of the largest medical practice groups in the country specializing in facial rejuvenation procedures.

Lifestyle Lift(C) chose the eCO2 after extensive research of laser platforms on the market. The eCO2's combination of safety, efficacy and value made it the best fit for the company's business model, its patients and its commitment to safety.

"We tested Lutronic's eCO2 on over one thousand patients and experienced the highest satisfaction among patients and physicians, as well as exceptional results," says Dr. David M. Kent, founder and CEO of Lifestyle Lift(C).

FDA-cleared since 2008 to treat wrinkles, furrows, fine lines, textural irregularities, pigmented lesions, vascular dyschromia and various other skin conditions for total skin resurfacing, the eCO2's ability to adjust depth penetration along with its proprietary Controlled Chaos scanning technology has resulted in a high level of efficacy with decreased discomfort, positioning it as a top player among multi-functional fractional CO2 laser platforms on the market - and making it an ideal choice for Lifestyle Lift(C).


The Lifestyle Lift(C) has become an increasingly popular alternative to traditional facelifts because it is safer, less expensive, and requires less recovery time than traditional facelifts. The Lifestyle Lift(C) procedure usually takes about one hour to perform, requires only local anesthesia, and can be done in the doctor's office in a relaxed atmosphere.

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
LASER09
ASLMS Conference

A Monumental Sharing of Laser Innovation

2009 ASLMS Annual Conference
April 1-5, 2009
Gaylord National™ Resort & Convention Center
near Washington D.C., April 1-5, 2009

ASLMS
SURGERY • TECHNOLOGY • MEDICINE™

14:26 - 15:14

Room3 Hakuho

SATURDAY OCTOBER 10

Laser

Chair : Shingo Wakamatsu [Japan]

O-131 Port-wine stains photodynamic therapy with intradermal administration of 5-aminolevulinic acid
Wei Li [China]

O-132 CO₂ fractional resurfacing with controlled chaos technique
Chang-Hun Huh [Korea]

O-133 New melasma treatment by collimated low fluence Q-switched Nd:YAG laser
Il-Hwan Kim [Korea]

O-134 Long-term evaluation of collagen and elastin following infrared (1100 to 1800 nm) irradiation
Yohei Tanaka [Japan]

O-135 Cosmetic effects and convenience of new 1550nm fractional Erbium-glass laser
Ga-Young Lee [Korea]

Won-serk, Kim

O-136 A multi-functional CO₂ fractional technology for comprehensive aesthetic treatments
Youjin Yang [Korea]

THE **AESTHETIC** GUIDE™
The Leading Aesthetic Practice Resource™

Beyond BOTOX
Home-Use Devices
Hand Rejuvenation

eCO2 Laser
Redefines
Fractional
Technology

For many years carbon dioxide (CO₂) lasers dominated laser resurfacing technology. Consumers were thrilled with the significant results, however, the associated risks and lengthy downtimes with this highly ablative treatment modality led to a decrease in interest. The patient's desire for less invasive treatments has influenced the skin rejuvenation field as a whole and skin resurfacing is no exception. Lutronic, Inc. (Princeton Junction, NJ) has combined the ablative properties of a CO₂ laser with fractional and scanning technology to create their eCO₂ devices. This unique combination provides a wide variety of treatment options and considerably reduces adverse effects.

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Global Standard
New Anti-aging Solution

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Effects of Ablative 10,600-nm Carbon Dioxide Fractional Laser Therapy on Suppurative Diseases of the Skin: A Case Series of 12 Patients

Sung Bin Cho, MD,¹ Jin Young Jung, MD,¹ Dong Jin Ryu, MD,¹ Sang Ju Lee, MD, PhD,² and Ju Hee Lee, MD, PhD^{1*}

¹Department of Dermatology and Cutaneous Biology Research Institute, Yonsei University College of Medicine, Seoul, Korea

²Yonsei Star Skin & Laser Clinic, Seoul, Korea

Background and Objectives: We have used an ablative 10,600-nm carbon dioxide fractional laser system (CO₂ FS) for suppurative diseases in order to attempt improvement. The purpose of our study was to demonstrate the effect of CO₂ FS on the course of inflammatory reactions in suppurative diseases.

Materials and Methods: We reviewed a total of 12 Korean patients with suppurative diseases of the skin who had a history of treatment failure with several therapeutic modalities as well as active and multiple inflammatory lesions at the time of CO₂ FS treatment.

Results: Improvement scores considering the number of suppurative lesions revealed that 3 of the 12 patients demonstrated clinical improvement of grade 4. Seven had clinical improvement of grade 3 and two showed improvement of grade 2. Improvement scores in severity were also evaluated; 2 of the 12 patients showed clinical improvement of grade 4. Six demonstrated clinical improvement of grade 3 and four had clinical improvement of grade 2. No patient showed a worsening of suppurative lesions.

Conclusion: Our observations demonstrated that the use of CO₂ FS did not make active suppurative lesions worse, and might have a therapeutic effect on suppurative diseases and their related scars. *Lasers Surg. Med.* 41:550–554, 2009. © 2009 Wiley-Liss, Inc.

Key words: fractional laser; carbon dioxide; acne vulgaris; folliculitis; furunculosis; hidradenitis suppurativa

INTRODUCTION

Fractional lasers have been proven to be effective and generally well tolerated for the treatment of various dermatologic diseases [1–3]. Although the non-ablative 1,550-nm erbium-doped fractional photothermolysis system (FPS) is the most popular fractional laser and yields decidedly reduced recovery times, the ablative 10,600-nm carbon dioxide fractional laser system (CO₂ FS) has recently been introduced and produces excellent clinical outcomes with fewer treatment sessions [4]. However, reports of the effects of fractional lasers usually consider acne scars rather than active lesions. Many dermatologists are also reluctant to use FPS or CO₂ FS for scars if patients present with active suppurative lesions.

A recent study described a patient with inflammatory acne vulgaris and acne scars who had successfully been treated with three sessions of combined 595-nm pulsed-dye laser with dynamic cooling device and FPS [5]. However, because the 595-nm pulsed-dye laser has proven therapeutic effects on acne vulgaris, the effect of FPS on active lesions could not be elucidated. In this report, we demonstrated the effect of CO₂ FS on the course of inflammatory reactions in suppurative diseases through a retrospective analysis of 12 patients.

MATERIALS AND METHODS

This retrospective study was performed in accordance with the ethical guidelines of the 1975 Declaration of Helsinki as reflected in approval by the Institutional Review Board of Severance Hospital, Yonsei University College of Medicine, Seoul, Korea. Twelve male patients (mean age: 22.3 years; range: 20–26 years; Fitzpatrick skin type: IV) who were treated using CO₂ FS for suppurative diseases of the skin with various clinical diagnoses, including inflammatory acne vulgaris ($n = 7$), chronic recurrent furunculosis ($n = 4$), and hidradenitis suppurativa (HS, $n = 1$), were reviewed in this study. Table 1 details subject characteristics. Patients with a history of treatment failure with several therapeutic modalities as well as active and multiple inflammatory lesions at the time of CO₂ FS treatment were enrolled. Subjects were excluded from this study if they had recently received systemic retinoids; CO₂ FS, FPS, 1,450-nm diode laser, 595-nm pulsed-dye laser, intense pulsed light, or photodynamic therapy within 6 months; or were treated with systemic and topical antibiotics, intradermal corticosteroid injection, incision and drainage, and surgical excision within 1 month.

Patients underwent 1–3 sessions (mean: 1.9 sessions) of CO₂ FS using a 10,600-nm Mosaic eCO₂™ (Lutronic Corporation, Goyang, Korea) laser. The Food and Drug Administration has approved the commercial release of

*Correspondence to: Ju Hee Lee, MD, PhD, Department of Dermatology and Cutaneous Biology Research Institute, Yonsei University College of Medicine, 250 Seongnam-ro, Seodamoon-gu, Seoul 120-752, Korea. E-mail: juhee@yuhs.ac


Accepted 30 June 2009
Published online in Wiley InterScience
(www.interscience.wiley.com).
DOI 10.1002/lsm.20802

Effects of Ablative 10,600-nm Carbon Dioxide Fractional Laser Therapy on Suppurative Diseases of the Skin: A Case Series of 12 Patients

Sung Bin Cho, MD,¹ Jin Young Jung, MD,¹ Dong Jin Ryu, MD,¹ Sang Ju Lee, MD, PhD,² and Ju Hee Lee, MD, PhD

¹Department of Dermatology and Cutaneous Biology Research Institute, Yonsei University College of Medicine, Seoul, Korea

²Yonsei Star Skin & Laser Clinic, Seoul, Korea

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Author Block Amina Husain, MD, Julie Woodward, MD.
Duke University Health System, Durham, NC, USA.

Author Disclosure Information: A. Husain, None.

Biography:

Dr. Husain is currently a fellow under Dr. Julie Woodward in oculoplastic surgery at Duke University Medical Center, Department of Ophthalmology

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Methods:

Retrospective chart review of patients undergoing bilateral lower lid laser resurfacing with the Lutronics eCO2 laser at the Duke Aesthetics Center.

- Charts were analyzed from January 2009 to July 2009.
- Photos and clinic visits were analyzed in terms of laser settings, Fitzpatrick grade, erythema, re-epithelization date, scarring, wrinkling, telangiectasias, dyschromia, pores, and patient satisfaction.
- Data was analyzed from day of procedure, post-operative week 1, month 1, month 3, and month 6.

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Ju Hee Lee, MD

*Department of Dermatology and Cutaneous Biology Research Institute,
Yonsei University College of Medicine, Seoul, Korea*

2009 IMCAS PARIS e-poster

IMCAS

International Master Course on Aging Skin
All the Leading Subjects on Aging Skin Treatments, Aesthetic Surgery, & More...

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1 pm - 2:15 pm POSTER SESSION

All posters will be displayed during the whole congress period at the Poster Area. A special Poster session is scheduled on Sunday. The author of the poster will have to stay next to his/her poster during this session in order to answer the participants questions.

POSTER NUMBER	PRACTITIONER	TOPIC
P805	Patrick BOGANOWICZ	Effects of hyaluronan on collagen synthesis and MMP-1 activity in stress induced premature senescent fibroblasts
P784	Sylvie BOISNIC	Evaluation of the tripolar technique using an ex-vivo model
P728	Francesca de ANGELIS	Non invasive full face skin rejuvenation based on combination of light sources and bipolar radiofrequency
P733	Karl Heinz EISELE	Is there a role for complexing proteins in pharmaceutical botulinum neurotoxin formulations?
P724	Fouad EL SAYED	Photopneumatic therapy in the treatment of recalcitrant acne and scars
P660	Nathalie FOURNIER	Cellulite treatment with photomology technology
P779	Nathalie FOURNIER	Efficient and painless laser skin tightening with pneumatic skin flattening
P714	Richard GENTILE	Technological innovation for facial rejuvenation procedures
P608	Julia GERSHKOVICH	A patient blinded prospective randomized controlled trial for correction of neolabial folds
P847	Mario GOISIS	Extracorporeal pulse-activation therapy: long-term effectiveness
P822	Zora HAVIAROVA	Body contouring – comparison of various ultrasound techniques
P665	Per HEDEN	Assessment of submuscular injection technique for breast enhancement
P668	Per HEDEN	Patient satisfaction after non-surgical breast enhancement
P723	Said HILTON	E-ART study result: an assessment of recovery time
P969	Klaus HOFFMANN	Volumizing effects of a smooth, highly cohesive, volumizing 20-mg/ml hyaluronic acid gel filler: prospective European study
P657	Oulkadi JAOUAD	Descriptive study about dual RF CONTRAGE in cellulite and abdominal flaccidity
P607	Joseph KABLIK	Evaluation of divinyl sulfone crosslinked acid dermal gel with lidocaine (DVS – NADGL) in a ginea pig skin sensitivity model
P717	Evangelos KERAMIDAS	A new classification and algorithm for treatment of breast ptosis with breast augmentation alone or simultaneous mastopexy and breast augmentation
P803	Franck LANDAT	Facial rejuvenation and contouring with a lasting polyacrylamide hydrogel
P877	Juan Pedro LA PUENTE	Resonant cavitation treatment of lipodystrophy results of the muscular research with 150 patients
P861	Hector LEAL SILVA	Male non invasive lipo-sculpturing using focused ultrasound measured by MRI
P852	Ju Hee LEE	Treatment of acne scars with 1550nm Er: Glass fractional laser and chemical reconstruction of skin scars (CROSS) method: a simultaneous split-face trial

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Jin Young Jung, MD et al.

Department of Dermatology and Cutaneous Biology Research Institute,
Department of Surgery and Department of Plastic and Reconstructive
Surgery, Yonsei University College of Medicine, Seoul, Korea

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- A 48-year-old woman with thyroidectomy scar (post-op 3 weeks).
(a) Before treatment. (b) 2 months after the 1st CO₂ FS treatment.
- Laser settings:
eCO2™ (50mJ/100 spots/cm²/static mode/1 pass/50% overlapping)

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- A 54-year-old woman with thyroidectomy scar (post-op 3 weeks).
(a) Before treatment. (b) 2 months after the 1st CO₂ FS treatment.
- Laser settings:
eCO2™ (50mJ/100 spots/cm²/static mode/1 pass/50% overlapping)

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- A 28-year-old woman with thyroidectomy scar (post-op 3 months).
(a) Before treatment. (b) 2 months after the 1st CO₂ FS treatment.
- Laser settings:
eCO₂TM (50mJ/100 spots/cm²/static mode/1 pass/50% overlapping)

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
Use of Intradermal Botulinum Toxin to Reduce Sebum Production and Facial Pore Size

Anil R. Shah MD

Clinical Instructor, Division of Facial Plastic Surgery, Department of Otolaryngology, University of Chicago, IL

ABSTRACT


Objective: Review the safety profile and subjective efficacy of intradermal botulinum toxin type A in facial pore size and sebum production.

Design: Retrospective analysis of 20 patients.

Results: Twenty consecutive patients with a single application of intradermal botulinum toxin type A were examined: Patients (17/20) noted an improvement in sebum production and a decrease in pores size at 1 month after injection. No complications were observed, and 17/20 patients were satisfied with the procedure.

Conclusion: Preliminary data suggests that intradermal botulinum toxin may play a role in decreasing sebum production. Further quantitative study may be necessary to determine effects of intradermal botulinum toxin on pore size.

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